



ALL IN ONE
SUPPLIES




CATALOGUE



Retractors
Mouth Gags



Wundhaken
Retractors
Divaricatori
Separadores
Écarteurs



JOSEPH
AI-BT-210/ AI-BT-214
160 mm




MINI-LANGENBECK
AI-BT-186/ AI-BT-189
160 mm

GILLIES
AI-BT-15
180 mm

GILLIES
AI-BT-16
180 mm

Wundhaken
Retractors
Divaricatori
Separadores
Écarteurs



1/1
Spitz
Sharp
Acuto
Agudo
Pointu

1/1
ADM-BT-38

1/1
ADM-BT-40

1/1
ADM-BT-41

1/1
ADM-BT-42

1/1
Stumpf
Blunt
Smusso
Romo
Mousse

1/1
ADM-BT-39

1/1
ADM-BT-44


1/1
ADM-BT-45

1/1
ADM-BT-46

ADM-BT-38/ADM-BT-46
160 mm

Feine wundhaken flexibel-spitz
Delicate retractors flexible-sharp
Divaricatori delicati flessibili-acuti
Separadores finos flexibles-agudos
Écarteurs fins pointu flexible

Feine wundhaken flexibel-stumpf
Delicate retractors flexible-blunt
Divaricatori delicati flessibili-smussi
Separadores finos flexibles-romos
Écarteurs fins mousse flexible



AI-BT-50
150 mm

AI-BT-51
150 mm

AI-BT-52
150 mm


AI-BT-54
150 mm

AI-BT-55
150 mm

AI-BT-56
150 mm



Wundhaken
Retractors
Divaricatori
Separadores
Écarteurs



AI-BT-70
10 mm

AI-BT-72
14 mm

AI-BT-74
15 mm

AI-BT-77
17 mm

DESMARRES
AI-BT-70/ AI-BT-77
140 mm



AI-BT-18

AI-BT-19


Spitz
Sharp
Acuto
Agudo
Pointu

Stumpf
Blunt
Smusso
Romo
Mousse

SEAN-MILLER
AI-BT-18/ AI-BT-19
160 mm




KOCHER
AI-BT-111/ AI-118
205 mm




AI-BT-117
19 x 12 mm

AI-BT-118
21 x 14 mm




HURD
AI-BT-1/ AI-BT-2
220 mm

Wundhaken
Retractors
Divaricatori
Separadores
Écarteurs







Wundhaken
Retractors
Divaricatori
Separadores
Écarteurs



AI-BT-265



AI-BT-264



DOCKHORN
AI-BT-134/ AI-BT-135
210 mm


AI-BT-134
40 x 12 mm



AI-BT-135
55 x 14 mm

AI-BT-264/ AI-BT-265
235 mm

Wundhaken
Retractors
Divaricatori
Separadores
Écarteurs





Wundhaken
Retractors
Divaricatori
Separadores
Écarteurs



OBWEGESER
AI-BT-152/ AI-BT-157
220 mm



OBWEGESER
AI-BT-92/ AI-BT-100
220 mm



AI-BT-152
25 x 7 mm

1/1




AI-BT-153
31 x 10 mm

1/1



AI-BT-154
42 x 11 mm

1/1



AI-BT-155
55 x 12 mm

1/1




AI-BT-156
70 x 14 mm

1/1




AI-BT-157
80 x 16 mm

1/1



AI-BT-92
25 x 7 mm

1/1



AI-BT-93
31 x 10 mm

1/1



AI-BT-94
42 x 11 mm

1/1



AI-BT-95
55 x 12 mm

1/1



AI-BT-96
70 x 14 mm

1/1




AI-BT-100
80 x 16 mm

1/1



Wundhaken
Retractors
Divaricatori
Separadores
Écarteurs



BOWLED-HENRY
AI-DO-175
170 mm




10 mm
AI-DO-204
170 mm


15 mm
AI-DO-204L
170 mm



AUSTIN
AI-DO-212
140 mm



1/1



AI-BT-33
28 x 10 mm


AI-BT-60
25 x 6 mm



AI-BT-34
30 x 14 mm



AI-BT-62
35 x 8 mm



AI-BT-35
33 x 14 mm



AI-BT-64
35 x 11 mm



AI-BT-36
40 x 10 mm



AI-BT-66
40 x 11 mm



AI-BT-43
40 x 13 mm



AI-BT-260
55 x 11 mm

LANGENBECK
AI-BT-33/ AI-BT-43
210 mm

KOCHER-LANGENBECK
AI-BT-60/ AI-BT-260
215 mm




Wundhaken
Retractors
Divaricatori
Separadores
Écarteurs


Rostfrei - Stainless Steel
Acciaio - Acero - Acer

Biegsam - Malleable
Malleabili - Maleables
Souples

ADM-OM-85
Kompletset mit 3 Spatel
Set complete with 3 spatulas
Set completo con 3 spatole
Set completos con 3 espátulas
Set complete avec 3 spatules
ADM-OM-80 + ADM-OM-81
ADM-OM-82




Wundhaken
Retractors
Divaricatori
Separadores
Écarteurs






Wundhaken
Retractors
Divaricatori
Separadores
Écarteurs



Mundsperrer
 Mouth gags
 Apribocca
 Abrebocas
 Écarteur buccal



Paar Einsätze aus Silikon für:
 Pair of silicone inserts for:
 Coppia di inserti in silicone per:
 Pareja de inserto de silicone para:
 Paire de inserts en silicon pour:

AI-OM-40 - AI-OM-41

Ø 8 mm


Ø 11 mm




ADM-OM-40R
 30 mm




AI-OM-41R
 30 mm



HEISTER
 AI-OM-62
 135 mm




ROSER KÖNIG
 AI-OM-60
 185 mm




Für Erwachsene - For adults
 Per adulti - Para adultos
 Pour adultes

Für Erwachsene - For adults
 Per adulti - Para adultos
 Pour adultes


Für kinder
 For children
 Per bambini
 Para niños
 Pour enfants



MOLT
 DOYEN-JANSEN
 AI-OM-41
 140 mm



MOLT
 DOYEN-COLLIN
 AI-OM-40
 120 mm




MOLT
 AI-OM-42
 100 mm


Mundsperrer
Mouth gags
Apribocca
Abrebocas
Écarteur buccal

Abhalter für gesamten Ober- bzw. Unterkiefer
Retractor for complete upper resp. Lower jaws
Retratore per arcata superiore e inferiore intera
Separadores para quijadas enteras inferiores y superiores
Écarteurs pour machoires complètes inférieures et supérieures


- AI-OM-74 Fig. 0 100 x 95 mm
- AI-OM-75 Fig. 1 115 x 105 mm
- AI-OM-76 Fig. 2 115 x 110 mm
- AI-OM-77 Fig. 3 125 x 115 mm



ORINGER
AI-OM-74/ AI-OM-77




MATY
AI-OM-97
70 mm




AI-OM-90
150 mm

Wundhaken
Retractors
Divaricatori
Separadores
Écarteurs






Wangenhalter
Cheek Retractors
Divaricatori per guance
Separador para mejillas
Écarteurs des joues




STERNBERG
AI-OM-50
160 mm



STERNBERG
AI-OM-52
160 mm




COLUMBIA
AI-OM-67
140 mm



Schwarz Verchromt
Black anodized
Cromato nero
Cromado en negro
Chromé noir

22 mm



BISHOP
AI DM-OM-61
130 mm



MINNESOTA
AI-OM-55
140 mm



CAWOOD-MINNESOTA
AI-OM-66
160 mm




MINNESOTA
AI-OM-55B
140 mm



CAWOOD-MINNESOTA
AI-OM-66B
160 mm

NEW
BLACK

Wundhaken
Retractors
Divaricatori
Separadores
Écarteurs




Ramus Retraktoren Recht
Ramus Retractor Right
Retrattore Ramus destro
Separadores Ramus derecho
Écarteur Ramus droites


Ramus Retraktoren Link
Ramus Retractor Left
Retrattore Ramus sinistro
Separadores Ramus izquierdo
Écarteur Ramus gauches

8 mm

10 mm



OBWEGESER
AI-DO-41
165 mm



OBWEGESER
AI-DO-142
165 mm




MISCH
AI-BT-158
210 mm




MISCH
AI-BT-159
210 mm

Mundsperrer
Mouth Gags
Apribocca
Abrebocas
Écarteur buccal



WHITEHEAD-JENNINGS

- AI-OM-136 90 mm
- AI-OM-137 110 mm
- AI-OM-138 130 mm
- AI-OM-139 150 mm




Komplett mit Rahmen, 3 Zungenspateln
und 2 beweglichen Seitenvälven
Complete with frame, 3 tongue depressors
and 2 adjustable lateral blades
Completo con telaio, 3 abbassalingua
e 2 valve laterali mobili
Completo con bastidor, 3 bajalenguas
y 2 valvas laterales ajustables
Complet avec cadre, 3 abaisse-langue
et 2 valves latérales réglables

DINGMANN
AI-OM-110
Set completo

DINGMANN


Zungenspatel
Tongue depressor
Abbassalingua
Bajalengua
Abaisse langue

Fig. 1




AI-OM-112
65 x 25 mm

Fig. 2




AI-OM-113
70 x 29 mm

Fig. 3



AI-OM-114
80 x 32 mm

Wundhaken
Retractors
Divaricatori
Separadores
Écarteurs




GREENE
AI-BT-131
225 mm

AI-BT-132
200 mm


WASSMUND
AI-BT-69
220 mm




Zungenspatel - Zungenzangen
Tongue Depressors - Tongue Holding Forceps
Abbassalingua - Pinze tiralingua
Bajalenguas - Pinzas tira-lengua
Abaisse langues - Pincés tire-langue




AI-OM-07
145 mm




KROENIG
AI-OM-16
150 mm




BUCHWALD
AI-OM-18
180 mm




BRUENINGS
AI-OM-04
190 mm





FRENZEL
AI-OM-19
165 mm



YOUNG
AI-OM-45
160 mm




COLLIN
AI-OM-10/ AI-OM-30




AI-OM-44

- 1 Paar Ersatzgummi
- 1 Pair of replacement rubber pads
- 1 Paio de gommini di ricambioi
- 1 Par de gomas de recambio
- 1 Paire de caoutchouc de rechange



AI-OM-30
165 mm



AI-OM-10
170 mm

Kette und Mundstutzen
 Chain and Mouth Props
 Catene e Apribocca
 Cadena y Abrebocas
 Chaîne et Ouvre Bouche



AI-KE-20
200 mm



AI-KE-40
400 mm

Dampfsterilisation bis 134°C
 Autoclavable to 134°C
 Autoclavabile a 134°C
 Autoclave hasta 134°C
 Autoclavable jusqu'à 134°C

Packung mit 2 Stück
 Package with 2 units
 Confezione da 2 pezzi
 Caja con 2 unidades
 Boîte avec 2 pièces

Schwarzer Gummi
 Black Rubber
 Gomma nera
 Caucho negro
 Caoutchouc noir



AI-OM-70
SMALL
Klein, für Kinder
Small, for children
Piccolo, per bambini
Pequeño, para niños
Petit, pour enfants



AI-OM-71
MEDIUM
Normal, für Erwachsene
Regular, for adults
Normale, per adulti
Regular, para adultos
Regulier, pour adultes



AI-OM-72
LARGE
Groß, für Erwachsene
Large, for adults
Grande, per adulti
Grande, para adultos
Grand, pour adultes



ALL IN ONE

SUPPLIES

Manufacturer & Exporter Of Surgical, Dental & Orthopedic Instruments.

+92-333-8740-629

info@aiosupplies.com | medisuraj@gmail.com

Bismillah Street, Kacha Shahab Pura, S.I.E. Sialkot Pakistan www.aiosupplies.com



ISO
9001:2008

